



Phenolic Resin Bedside Table

AMIGO 1

We Provide Customized Overall Solutions
for Medical Equipment

PRODUCT SPECIFICATIONS

Folded size (LxWxH)	475x475x800 mm
---------------------	----------------

Drawer	1pc
Cabinet	1pc
Casters	4pcs

QUALITY CONTROL

Handle

The handle is made of aluminum alloy with surface anodized, which is firm, rust-proof and easy to clean.



Casters

Casters is using medical silent material, double side casters, stable and noiseless. 2 casters with brake, lockable at 360 degree, 2 casters free and swivel to any direction.




Frame

Aluminum alloy 4 column frame, stable and firm.



Material

Table body 3 side walls are made of powder coated steel, table top, drawer and door are made of phenolic compact board with rounding treatment, anti-acid, anti-alkali, waterproof and anti-seepage, easy to clean and flame retardant for 4 hours.



Structure

A drawer, a cabinet, a storage open case in the middle.




Sliding Track

The sliding track adopts silent damping track, which can effectively prevent the drawer from sliding out in the course of moving the cart.



Color

Option: color (color samples) .



Comprehensive Hospital Design Solutions




Toruscare has a complete industry chain targeting the global market.



Toruscare provides overall healthcare solutions.



Toruscare provides services for full coverage of target audiences.



Design and Development Excellence

At Toruscare, our strength lies in our robust design and development capabilities. Our dedicated team of engineers and designers continually push the boundaries of innovation, ensuring that our products are at the forefront of medical technology. From conceptualization to prototyping, we meticulously craft each product to meet the highest standards of functionality, durability, and ease of use.



Uncompromising Quality Control

Toruscare is deeply committed to quality, implementing stringent control measures at every manufacturing stage. Our state-of-the-art facilities, equipped with advanced testing and inspection tools, ensure each product meets rigorous standards. This dedication to excellence has earned the trust and loyalty of healthcare providers worldwide.



Global Reach, Local Touch

As a global company, Toruscare understands the diverse needs of different markets. We pride ourselves on our ability to deliver customized solutions that cater to the unique requirements of healthcare providers across the world. Our extensive distribution network and strategic partnerships ensure timely delivery and exceptional after-sales support, making us a trusted partner in the global healthcare community.



Care for All

TORUSCARE

Toruscare Sp. z o.o

Adress: UL. SPACEROWA 1C, 41-253 CZELADŹ, POLAND

Email: export@toruscare.com

Phone: +351 910200400