



Electric Examination Table

# GANAS 1

We Provide Customized Overall Solutions  
for Medical Equipment



# PRODUCT SPECIFICATIONS

External Size (LxWxH)	1810x680x(480-850) mm
Back-Rest Adjustment	0-50°
Weighh	150 kg

Pillow	
Motor	
Mechanical Lever to Adjust Back Rest	
2" Casters	
Hand Controller	
Paper Roll Holder	

1pc
1pc
1pc
2pcs
1pc
1pc

# POSTURE FUCTION

Back-Rest



Height Adjustment



# HOSPITAL-GRADE EXAMINATION BED

Flame retardant, yellow resistance, mildew resistance, bending resistance, 100,000 times wear resistance.



1. 2 pcs 2' Hidden castors, 2 pcs 2' universal wheels.
2. Back adjust can choose ratchet or motor.
3. Stainless steel foot bar system, without wear or tear after 5000 times use.
4. Flame retardant, 8 cm thickness, high-density foam, density 30 kg/m<sup>3</sup>.
5. Engineering ABS material plug, foot pad.
6. 2 Colors as option, black and blue-white.



## Accessories

### Leather

Flame retardant, yellow resistance, mildew resistance, bending resistance, 100,000 times wear resistance.



### Sponge

Flame retardant, 10 cm thickness, high-density foam.



## Painting

11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



## Back adjust

Back adjust can choose ratchet or motor.



## Castor

2 pcs 2-inch universal castors, with foot pedal.



## Castor

2 pcs 2-inch Hidden castors.



## Foot pad

Engineering ABS material plug, foot pad.



## Hand controller

Hand controller, keep qualified after 2000 times use.





# Comprehensive Hospital Design Solutions




Toruscare has a complete industry chain targeting the global market.



Toruscare provides overall healthcare solutions.



Toruscare provides services for full coverage of target audiences.



## Design and Development Excellence

At Toruscare, our strength lies in our robust design and development capabilities. Our dedicated team of engineers and designers continually push the boundaries of innovation, ensuring that our products are at the forefront of medical technology. From conceptualization to prototyping, we meticulously craft each product to meet the highest standards of functionality, durability, and ease of use.



## Uncompromising Quality Control

Toruscare is deeply committed to quality, implementing stringent control measures at every manufacturing stage. Our state-of-the-art facilities, equipped with advanced testing and inspection tools, ensure each product meets rigorous standards. This dedication to excellence has earned the trust and loyalty of healthcare providers worldwide.



## Global Reach, Local Touch

As a global company, Toruscare understands the diverse needs of different markets. We pride ourselves on our ability to deliver customized solutions that cater to the unique requirements of healthcare providers across the world. Our extensive distribution network and strategic partnerships ensure timely delivery and exceptional after-sales support, making us a trusted partner in the global healthcare community.



Care for All

**TORUSCARE**

Toruscare Sp. z o.o

Adress: UL. SPACEROWA 1C, 41-253 CZELADŹ, POLAND

Email: [export@toruscare.com](mailto:export@toruscare.com)

Phone: +351 910200400