



Electric Examination Table

## GANAS 3

We Provide Customized Overall Solutions  
for Medical Equipment



## PRODUCT SPECIFICATIONS

External Size (LxW)	1960x680 mm	Rated voltage and frequency	AC220V, 50Hz
Height Adjustment	470-910 mm	Input power	500 VA
Head up/down	88°/40°		
Leg up/down	50°/45°		
Net weight	68 kg		
Safe working load	160 kg		

Adjustable arm support	1pc
Motor	1pc
Gas spring	2pcs
3" Castor	4pcs
Brake Pedal	4pcs
Hand Controller	1pc
Manual Operator	1pc



# KEY FEATURES

Features a high-quality motor, overall iron frame structure, smooth and quiet height and angle adjustments.

Equipped with premium-grade caster lift switches for easy transition between fixed and mobile modes; smooth and quiet movement.



High-quality steel frame for sturdiness and stability, ensuring safe use.

Antibacterial leather and high-density slow-rebound sponge provide comfort while meeting antibacterial standards.



# Accessories

## Framework

Bed frame surface painted with 1.5mm powder coated steel, seamless welded molding by robot, strong corrosion resistance, environmentally-friendly, and bright surface maintain an overall beautiful appearance. With this structure, the whole bed maintains more stable and makes patients feel more safe.



## Mattress

The bed surface is made of the integrated high-end composite leather and hand-made cotton mattress, of which the curve design is ergonomic, elastic and soft, breathable, environmentally-friendly, easy to be cleaned, anti-bacterial, impact-resistant.



## Control system

Manual air levers with two air springs to adjust the back and leg plate angle provide users with a quiet and safe environment during diagnosis and avoid interference.



## Braking system

Equipped with Pedal-controlled 3 castors: 4pcs, brake pedals in 2 functions-locked and universal: 4 pcs, convenient to move and fix. 360° spinning wheels polyurethane PU-strong, durable, non-slip, strong in pressure resistance, and good in stability. Move freely when sitting on the bed, saving energy and effort.



## Hand controller

Equipped with hook-type hand controller, easy to press and adjust optimal height. The hook on the back is convenient to take down for patients and nurses.



## Matching Pillow

Neoteric design, narrow on the top and wide on the bottom, can be disassembled as required, easy to use. Its well-shaped and elegant design matches with the modern aesthetic, making it harmoniously integrates into the environment.



## Round Covered Edge

The bed platform was covered in 360°. The soft-touch of surface would be safer for users.



## Design

The bed surface with two folds of three sections ergonomic design, can realize elevation adjustment according to the patient's height, achieving personalized design to help diagnosis. Bed and chair dual.use, high load-bearing, durable and stable, fearless of heavy pressure impact, to ensure the safety of use.



## Optional Foot-pedal Controller

Available to operate the entire bed, which can be adjusted to the appropriate height by lightly stepping on it.



# Comprehensive Hospital Design Solutions




Toruscare has a complete industry chain targeting the global market.



Toruscare provides overall healthcare solutions.



Toruscare provides services for full coverage of target audiences.



## Design and Development Excellence

At Toruscare, our strength lies in our robust design and development capabilities. Our dedicated team of engineers and designers continually push the boundaries of innovation, ensuring that our products are at the forefront of medical technology. From conceptualization to prototyping, we meticulously craft each product to meet the highest standards of functionality, durability, and ease of use.



## Uncompromising Quality Control

Toruscare is deeply committed to quality, implementing stringent control measures at every manufacturing stage. Our state-of-the-art facilities, equipped with advanced testing and inspection tools, ensure each product meets rigorous standards. This dedication to excellence has earned the trust and loyalty of healthcare providers worldwide.



## Global Reach, Local Touch

As a global company, Toruscare understands the diverse needs of different markets. We pride ourselves on our ability to deliver customized solutions that cater to the unique requirements of healthcare providers across the world. Our extensive distribution network and strategic partnerships ensure timely delivery and exceptional after-sales support, making us a trusted partner in the global healthcare community.



Care for All

**TORUSCARE**

Toruscare Sp. z o.o

Adress: UL. SPACEROWA 1C, 41-253 CZELADŹ, POLAND

Email: [export@toruscare.com](mailto:export@toruscare.com)

Phone: +351 910200400