



New All-In-One Computer Cart

IMPETU LT1

We Provide Customized Overall Solutions
for Medical Equipment

PRODUCT SPECIFICATIONS

Size (LxWxH)	7490×460×780-1030 mm
--------------	----------------------

Big Basket	1pc
Tray	1pc
Castor	4pcs

QUALITY CONTROL

Column

Aluminum alloy column.



Top surface

ABS material top, easy to clean and durable to scratch. Cover is free to open, laptop can be put in easily.



Compartment

Compartment can storage any laptop size under 22", load capacity: 10kg. LED keyboard lights illuminate the keyboard.



Mouse holder

Double sided telescopic holder for mouse, left or right hand operating are both available.



Base

Base is made of ABS material, featured with both stability and easy-maintenance, with foot pedal to adjust height.



Castor

Castors are made of medical material, double sided, stable and noiseless. 2 front casters with brake, lockable at 360 degree, 2 rear casters are flexible and directional.



Computer holder

With basket on the back, easy to place goods, Environmentally plastic is non-toxic and easy to clean.



Battery

Optional configuration with high-performance lithium battery (with LED Keyboard light), constantly supply stable power source, safe and convenient.



Comprehensive Hospital Design Solutions




Toruscare has a complete industry chain targeting the global market.



Toruscare provides overall healthcare solutions.



Toruscare provides services for full coverage of target audiences.



Design and Development Excellence

At Toruscare, our strength lies in our robust design and development capabilities. Our dedicated team of engineers and designers continually push the boundaries of innovation, ensuring that our products are at the forefront of medical technology. From conceptualization to prototyping, we meticulously craft each product to meet the highest standards of functionality, durability, and ease of use.



Uncompromising Quality Control

Toruscare is deeply committed to quality, implementing stringent control measures at every manufacturing stage. Our state-of-the-art facilities, equipped with advanced testing and inspection tools, ensure each product meets rigorous standards. This dedication to excellence has earned the trust and loyalty of healthcare providers worldwide.



Global Reach, Local Touch

As a global company, Toruscare understands the diverse needs of different markets. We pride ourselves on our ability to deliver customized solutions that cater to the unique requirements of healthcare providers across the world. Our extensive distribution network and strategic partnerships ensure timely delivery and exceptional after-sales support, making us a trusted partner in the global healthcare community.



Care for All

TORUSCARE

Toruscare Sp. z o.o

Adress: UL. SPACEROWA 1C, 41-253 CZELADŹ, POLAND

Email: export@toruscare.com

Phone: +351 910200400