



Overbed Table

MESITA 1

We Provide Customized Overall Solutions
for Medical Equipment

PRODUCT SPECIFICATIONS

External size (LxWxH)	800×395×870-1170mm	Air Spring for Height Adjustment	1pc
		Phenolic Resin Table Top	1pc

QUALITY CONTROL

Material

Powder coated steel base, powder coated steel column.



Table Top

Table top is made of phenolic resin laminate with PVC cover.



Height Adjustment

Height adjustment by gas spring, used single-handed according to the convenience of the user.



Casters

Easy to install, more stable and firmer

High quality medical caster, 1.5 inch /2 inch mute wheel.



Matching Design

It can be used with the bedside table, which is easy to disassemble.
The combination of dining table and cabinet functions can serve two purposes and save space.




Comprehensive Hospital Design Solutions




Toruscare has a complete industry chain targeting the global market.



Toruscare provides overall healthcare solutions.



Toruscare provides services for full coverage of target audiences.



Design and Development Excellence

At Toruscare, our strength lies in our robust design and development capabilities. Our dedicated team of engineers and designers continually push the boundaries of innovation, ensuring that our products are at the forefront of medical technology. From conceptualization to prototyping, we meticulously craft each product to meet the highest standards of functionality, durability, and ease of use.



Uncompromising Quality Control

Toruscare is deeply committed to quality, implementing stringent control measures at every manufacturing stage. Our state-of-the-art facilities, equipped with advanced testing and inspection tools, ensure each product meets rigorous standards. This dedication to excellence has earned the trust and loyalty of healthcare providers worldwide.



Global Reach, Local Touch

As a global company, Toruscare understands the diverse needs of different markets. We pride ourselves on our ability to deliver customized solutions that cater to the unique requirements of healthcare providers across the world. Our extensive distribution network and strategic partnerships ensure timely delivery and exceptional after-sales support, making us a trusted partner in the global healthcare community.



Care for All

TORUSCARE

Toruscare Sp. z o.o

Adress: UL. SPACEROWA 1C, 41-253 CZELADŹ, POLAND

Email: export@toruscare.com

Phone: +351 910200400