



Electric Homecare Bed

HOGA 1

We Provide Customized Overall Solutions
for Medical Equipment

 www.toruscare.com

TORUSCARE

Care for All



PRODUCT SPECIFICATIONS

External Size	2140x1055 mm	Lowest & Highest Position	190-610 mm
Mattress platform	2000x900 mm	Safe working load	215 kg
Back-rest Adjustment	0-72°	Mattress size	L2000×900 mm/800 mm
Leg-rest Adjustment	0-28°		
Reverse Trendelenburg	0-12°		
Trendelenburg	0-12°		

Linear Motor	4pcs
Wooden side rails	1set
Head & Foot Board	1set
Castors	4pcs
Remote controller	1pc
Emergency powersystem	Optional
Trapeze Pole	Optional




GREAT MULTI-CARE

Good night, better days

SUPER LOW POSITION

Lowest height 190mm from the ground, to prevent accidental falls from the bed.



KEY FEATURES

1. Wide range of height adjustment (lowest position-19 cm);
2. 4-column construction;
3. Easy and safe to fold and transport;
4. Lowered position of side rails;
5. Various wood colours to choose from;
6. Maximum working load-215 kg;
7. 2-year warranty;

Adjust design to individual needs

Wide range of wood colours:



Beech



Light oak



Maple

Dark walnut

Gold oak


Sonoma Tabac

Palermo oak light

Accessories

Transport Stand

The whole bed can be folded for transportation. Safe and comfortable to transport. Easy to place and can be folded when not in use.



Head & Foot Board

The head and foot boards are made of laminate and solid beech wood, providing strong stability while being lightweight.



Siderail

Wooden-edged guardrails maximize patient safety and prevent falls. For comfort, the guardrail is lowered to facilitate patients to stand up and sit on the bed.



Remote Controller

The 10-button handset is equipped with a locking function to prevent accidental touches.



Comprehensive Hospital Design Solutions




Toruscare has a complete industry chain targeting the global market.



Toruscare provides overall healthcare solutions.



Toruscare provides services for full coverage of target audiences.



Design and Development Excellence

At Toruscare, our strength lies in our robust design and development capabilities. Our dedicated team of engineers and designers continually push the boundaries of innovation, ensuring that our products are at the forefront of medical technology. From conceptualization to prototyping, we meticulously craft each product to meet the highest standards of functionality, durability, and ease of use.



Uncompromising Quality Control

Toruscare is deeply committed to quality, implementing stringent control measures at every manufacturing stage. Our state-of-the-art facilities, equipped with advanced testing and inspection tools, ensure each product meets rigorous standards. This dedication to excellence has earned the trust and loyalty of healthcare providers worldwide.



Global Reach, Local Touch

As a global company, Toruscare understands the diverse needs of different markets. We pride ourselves on our ability to deliver customized solutions that cater to the unique requirements of healthcare providers across the world. Our extensive distribution network and strategic partnerships ensure timely delivery and exceptional after-sales support, making us a trusted partner in the global healthcare community.



Care for All

TORUSCARE

Toruscare Sp. z o.o

Adress: UL. SPACEROWA 1C, 41-253 CZELADŹ, POLAND

Email: export@toruscare.com

Phone: +351 910200400