



Super Low Children Bed

HOGA CHD 1

We Provide Customized Overall Solutions
for Medical Equipment



PRODUCT SPECIFICATIONS

Bed dimensions	2040x926 mm	Safe working load	215 kg
Bed platform	1900x780 mm	Height adjustment	268-688 mm
Back rest	0-72°		
Leg rest	0-28°		
Reverse Trendelenburg	0-12°		
Mattress size	1900x780 mm		

Linear Motor	4pcs
Wooden side rails	1set
Head & Foot Board	1set
Castors	4pcs
Remote controller	1pc
Trapeze Pole	Optional



GREAT MULTI-CARE



268 mm

SUPER LOW POSITION

Lowest height 268mm from the ground, preventing children from accidentally falling off the bed.

KEY FEATURES

1. Wide range of height adjustment (lowest position-268mm);
2. Bright and eye-catching colours;
3. 4-section bed;
4. Easy and safe to fold and transport;
5. 3 side rails for maximum safety;
6. Handset grip located out of patient's reach;
7. Optional accessories;



Adjust design to individual needs

Wide range of wood colours:



Beech



Light oak



Accessories

Transport Stand

The whole bed can be folded for transportation. Safe and comfortable to transport. Easy to place and can be folded when not in use.



Head & Foot Board

The head and foot boards are made of laminate and solid beech wood, providing strong stability while being lightweight.



Siderail

Wooden-edged guardrails maximize patient safety and prevent falls. For comfort, the guardrail is lowered to facilitate patients to stand up and sit on the bed.



Remote Controller

The 10-button remote controller has a lock function to prevent accidental touches.



Comprehensive Hospital Design Solutions




Toruscare has a complete industry chain targeting the global market.



Toruscare provides overall healthcare solutions.



Toruscare provides services for full coverage of target audiences.



Design and Development Excellence

At Toruscare, our strength lies in our robust design and development capabilities. Our dedicated team of engineers and designers continually push the boundaries of innovation, ensuring that our products are at the forefront of medical technology. From conceptualization to prototyping, we meticulously craft each product to meet the highest standards of functionality, durability, and ease of use.



Uncompromising Quality Control

Toruscare is deeply committed to quality, implementing stringent control measures at every manufacturing stage. Our state-of-the-art facilities, equipped with advanced testing and inspection tools, ensure each product meets rigorous standards. This dedication to excellence has earned the trust and loyalty of healthcare providers worldwide.



Global Reach, Local Touch

As a global company, Toruscare understands the diverse needs of different markets. We pride ourselves on our ability to deliver customized solutions that cater to the unique requirements of healthcare providers across the world. Our extensive distribution network and strategic partnerships ensure timely delivery and exceptional after-sales support, making us a trusted partner in the global healthcare community.



Care for All

TORUSCARE

Toruscare Sp. z o.o

Adress: UL. SPACEROWA 1C, 41-253 CZELADŹ, POLAND

Email: export@toruscare.com

Phone: +351 910200400