



Rehabilitation Bed

## HOGA XL

We Provide Customized Overall Solutions  
for Medical Equipment



## PRODUCT SPECIFICATIONS

External Size	2195x1345 mm	Lowest & Highest Position	230-660 mm
Mattress platform	2010x1200 mm	Safe working load	320 kg
Back-rest Adjustment	0-70°	Mattress size	L2000 × 1200 mm
Leg-rest Adjustment	0-24°		
Reverse Trendelenburg	0-12°		
Trendelenburg	0-12°		

Linear Motor	4pcs
Wooden side rails	1set
Head & Foot Board	1set
Castors	4pcs
Remote controller	1pc
Emergency power system	Optional
Trapeze Pole	Optional



# GREAT MULTI-CARE

Good night, better days

## MAXIMUM PATIENT WEIGHT

The extra stable and sturdy frame allows a safe load of up to 320 kg.



## KEY FEATURES

1. Reinforced structure increases load capacity;
2. Bed width increased to 120cm;
3. Bed position lowered;
4. Maximum working load 372kg;
5. Various colors available;
6. Control bed functions with an ergonomic remote control;
7. Possibility of locking function;
8. 2-year warranty;

## Adjust design to individual needs

Wide range of wood colours:



Beech



Light oak





# Accessories

## Transport Stand

The whole bed can be folded for transportation. Safe and comfortable to transport. Easy to place and can be folded when not in use.



## Head & Foot Board

The head and foot boards are made of laminate and solid beech wood, providing strong stability while being lightweight.



## Siderail

Full-length heavy wooden guardrails (adjustable in height) come standard with this bed, providing protection up to 370 mm above the mattress platform.



## Remote Controller

The 10-button handset is equipped with a locking function to prevent accidental touches.




# Comprehensive Hospital Design Solutions




Toruscare has a complete industry chain targeting the global market.



Toruscare provides overall healthcare solutions.



Toruscare provides services for full coverage of target audiences.



## Design and Development Excellence

At Toruscare, our strength lies in our robust design and development capabilities. Our dedicated team of engineers and designers continually push the boundaries of innovation, ensuring that our products are at the forefront of medical technology. From conceptualization to prototyping, we meticulously craft each product to meet the highest standards of functionality, durability, and ease of use.



## Uncompromising Quality Control

Toruscare is deeply committed to quality, implementing stringent control measures at every manufacturing stage. Our state-of-the-art facilities, equipped with advanced testing and inspection tools, ensure each product meets rigorous standards. This dedication to excellence has earned the trust and loyalty of healthcare providers worldwide.



## Global Reach, Local Touch

As a global company, Toruscare understands the diverse needs of different markets. We pride ourselves on our ability to deliver customized solutions that cater to the unique requirements of healthcare providers across the world. Our extensive distribution network and strategic partnerships ensure timely delivery and exceptional after-sales support, making us a trusted partner in the global healthcare community.



Care for All

**TORUSCARE**

Toruscare Sp. z o.o

Adress: UL. SPACEROWA 1C, 41-253 CZELADŹ, POLAND

Email: [export@toruscare.com](mailto:export@toruscare.com)

Phone: +351 910200400