



Rehabilitation Therapy  
Emergency Trolley  
Medical Carts  
Pre Hospital Rescue  
Supplies

Emergency Trolley

## IMPETU EM1

We Provide Customized Overall Solutions  
for Medical Equipment

# PRODUCT SPECIFICATIONS

Overall Dimensions	790W X 610D X 1000H mm
Defibrillator Stand	380×340mm
Work Board	340×250mm
Cpr Board	475×595×8.5mm
Packing Size	820×650×1150 mm

Product Weight	0-35° (±10°)
Caster	5"Two Casters With Brake
Work Board Size	340×250mm
Bracket Size	380×340mm

20L Trash Can With Lid
Oxygen Bottle Rack
Adjustable Iv Stand
Medicine Grid Horizontal And Vertical
Slide Rail

## FEATURES AND HIGHLIGHTS

The ABS material is an innovative engineering plastic.Antibacterial, moisture-proof, anti-corrosion, wear-resistant, non-toxic and tasteless, never rust. Washable, easy cleaning.

- The overall design and production are exquisite, the styling and color matching are reasonable, installation is simple, and the structure is stable and durable.
- Baosteel from Fortune Global 500.Panasonic Robotic ensure 360° full smooth welding.
- There are transparent acrylic medicine grids inside the drawer, which can be used to classify and place medicines, and the medicine grids can be freely combined.



# QUALITY CONTROL

## IV stand

Adjustable aluminum alloy IV stand.



## Caster

At the bottom, it uses two 5-inch silent brake, anti-winding, high-end, silent casters, allowing for flexible movement and a high load-bearing capacity.



## Drawer slide

The drawer uses high-strength three-section silent damping slides. When the drawer is closed to the final distance, hydraulic pressure is used to slow down the speed and reduce the impact force.



## Working table

The main tabletop can have a replaceable soft pad that is dust-proof and easy to clean. The vehicle is equipped with a hidden sliding writing board that can be used for writing.



## Sealed handle

The trolley is equipped with a sealed handle, which can lock five drawers at the same time.



## Oxygen bottle rack

ABS oxygen bottle bracket on the left side.



## Trash can

Clamshell ABS trash can on the right side,

Specification: 20L



## CPR board

The ABS CPR board is installed on the back of the trolley body, which can be used during defibrillation rescue.

CPR board size: 475 \* 595 \* 8.5mm



## POSTURE FUNCTION





# Comprehensive Hospital Design Solutions




Toruscare has a complete industry chain targeting the global market.



Toruscare provides overall healthcare solutions.



Toruscare provides services for full coverage of target audiences.



## Design and Development Excellence

At Toruscare, our strength lies in our robust design and development capabilities. Our dedicated team of engineers and designers continually push the boundaries of innovation, ensuring that our products are at the forefront of medical technology. From conceptualization to prototyping, we meticulously craft each product to meet the highest standards of functionality, durability, and ease of use.



## Uncompromising Quality Control

Toruscare is deeply committed to quality, implementing stringent control measures at every manufacturing stage. Our state-of-the-art facilities, equipped with advanced testing and inspection tools, ensure each product meets rigorous standards. This dedication to excellence has earned the trust and loyalty of healthcare providers worldwide.



## Global Reach, Local Touch

As a global company, Toruscare understands the diverse needs of different markets. We pride ourselves on our ability to deliver customized solutions that cater to the unique requirements of healthcare providers across the world. Our extensive distribution network and strategic partnerships ensure timely delivery and exceptional after-sales support, making us a trusted partner in the global healthcare community.



Care for All

**TORUSCARE**

Toruscare Sp. z o.o

Adress: UL. SPACEROWA 1C, 41-253 CZELADŹ, POLAND

Email: [export@toruscare.com](mailto:export@toruscare.com)

Phone: +351 910200400