

# Electric Obstetric Bed



Electric Obstetric Bed


## ROSAWIDE 2

We Provide Customized Overall Solutions  
for Medical Equipment

 [www.toruscare.com](http://www.toruscare.com)

**TORUSCARE**

Care for All



## PRODUCT SPECIFICATIONS

Exteral Size (LxW)	2050x960 mm(±30mm)
Height Adjustment	640-1000 mm(±20mm)
Back-Rest Adjustment	0-70°
Reverse Trendelenburg	0-10° (±2°)
Trendelenburg	0-10° (±2°)
Leg-rest Adjustment	-5°-65°, leg section is removable

Weight capacity	200 KG
-----------------	--------

Linak Motor
5" Double Side Castors, Central Locking
Central Locking Pedal
ABS Head Board With Safe Lock
ABS Side Rail With Angle Indicator
Hand Grip
Leg Holder
Basin
Hand Controller
Power Line
Leg strip
IV Pole Provision

3pcs
4pcs
2pcs
1pc
1set
1pair
1pair
1pc
1pc
1pc
1pair
2pcs




## Material

The base seat and height adjustment system are made of powder coated steel, high gloss, anti corrosion, easy to clean.



## Bed surface

Bed surface by plastic steel, made by virgin material. The cover is detachable, and easy to clean.



## Head Board and Guardrail

The Head board and side rail made by virgin ABS material. The head board is movable.



## Hand controller

Equipped with one hand controller.



## Handle

The hand bar can be dismantled and movable.



## Lithotomy rod

The legsupport can be dismantled and movable.



## Caster

5 inch central controlled caster for moving obstetric table.  
TPR tire, quiet and wear-resistant, can be washed and sterilized  
by high temperature, save anti-winding hard shell, united forming  
without bolts.



## Accessory

Equipped with stainless steel sewage basin and other various  
accessories.



## Mattress

The mattress adopts high density memory sponge and is once  
formed.  
Mattress is movable.



## Foot section

Foot section can be retracted, convenient for gynecological and  
delivery surgical.



# Comprehensive Hospital Design Solutions




Toruscare has a complete industry chain targeting the global market.



Toruscare provides overall healthcare solutions.



Toruscare provides services for full coverage of target audiences.



## Design and Development Excellence

At Toruscare, our strength lies in our robust design and development capabilities. Our dedicated team of engineers and designers continually push the boundaries of innovation, ensuring that our products are at the forefront of medical technology. From conceptualization to prototyping, we meticulously craft each product to meet the highest standards of functionality, durability, and ease of use.



## Uncompromising Quality Control

Toruscare is deeply committed to quality, implementing stringent control measures at every manufacturing stage. Our state-of-the-art facilities, equipped with advanced testing and inspection tools, ensure each product meets rigorous standards. This dedication to excellence has earned the trust and loyalty of healthcare providers worldwide.



## Global Reach, Local Touch

As a global company, Toruscare understands the diverse needs of different markets. We pride ourselves on our ability to deliver customized solutions that cater to the unique requirements of healthcare providers across the world. Our extensive distribution network and strategic partnerships ensure timely delivery and exceptional after-sales support, making us a trusted partner in the global healthcare community.





Care for All

**TORUSCARE**

Toruscare Sp. z o.o

Adress: UL. SPACEROWA 1C, 41-253 CZELADŹ, POLAND

Email: [export@toruscare.com](mailto:export@toruscare.com)

Phone: +351 910200400