




IV Stand

PERO 1

We Provide Customized Overall Solutions
for Medical Equipment

PRODUCT SPECIFICATIONS

Height adjustment	1.2-2m
Weight	6.3kg
Column diameter	SS Outer pipe $\Phi 32\text{mm}$, inner pipe $\Phi 19\text{mm}$
Storage box size	$\phi 17\text{cm}$
Storage box load capacity	5kg
Base diameter	500mm



QUALITY CONTROL

Design

The appearance design is simple and elegant, lightweight, and easy to move.



Material

The main body is made of high-quality 201 stainless steel tubes, with good load-bearing capacity and not easily prone to rust.



Base

The steel plate base is equipped with an ABS dust cover, easy to clean, slip-resistant, and won't tip over.




Wheels

Five 3-inch wheels, flexible mobility, 2 wheels equipped with brakes.



Hooks

Four stainless steel hooks on the top that can be used by multiple individuals.



Name card

With a department name card for easy identification, which can also be used to attach medication documents.



Storage box

The ABS round storage box with a pull ring design allows for easy placement and convenient movement.



Plastic hooks

Add adjustable height plastic hooks at the bottom to hang other items besides medicines, making it more convenient for medical staff and patients.




ABS dust cover

Equipped with an ABS dust cover to prevent dust accumulation, ensuring a cleaner and more hygienic environment.



Height adjustment

Plastic knobs are used for height adjustment, allowing for free and flexible height customization. Add a height scale for easy height marking.



Comprehensive Hospital Design Solutions




Toruscare has a complete industry chain targeting the global market.



Toruscare provides overall healthcare solutions.



Toruscare provides services for full coverage of target audiences.



Design and Development Excellence

At Toruscare, our strength lies in our robust design and development capabilities. Our dedicated team of engineers and designers continually push the boundaries of innovation, ensuring that our products are at the forefront of medical technology. From conceptualization to prototyping, we meticulously craft each product to meet the highest standards of functionality, durability, and ease of use.



Uncompromising Quality Control

Toruscare is deeply committed to quality, implementing stringent control measures at every manufacturing stage. Our state-of-the-art facilities, equipped with advanced testing and inspection tools, ensure each product meets rigorous standards. This dedication to excellence has earned the trust and loyalty of healthcare providers worldwide.



Global Reach, Local Touch

As a global company, Toruscare understands the diverse needs of different markets. We pride ourselves on our ability to deliver customized solutions that cater to the unique requirements of healthcare providers across the world. Our extensive distribution network and strategic partnerships ensure timely delivery and exceptional after-sales support, making us a trusted partner in the global healthcare community.



Care for All

TORUSCARE

Toruscare Sp. z o.o

Adress: UL. SPACEROWA 1C, 41-253 CZELADŹ, POLAND

Email: export@toruscare.com

Phone: +351 910200400