




IV Stand

PERO 2

We Provide Customized Overall Solutions
for Medical Equipment

PRODUCT SPECIFICATIONS


Size (LxWxH)	500×490×(1310-2170)mm
Height adjustment	1310-2170mm
Inner tube	Φ16mm stainless steel round pipe
Outer tube	Φ25mm stainless steel round pipe
Base	20×40×2mm carbon steel tube



QUALITY CONTROL

Structure

Stainless steel telescopic with 4 hooks, plastic knob made of ABS. spring inside the tube.H Shaped feet, made of powder coated square tube.



Loading weight

Loading weight for each hook:2kg



Castor

Powder coated feet with 4 Castors, double side, two with brake. diameter:2inch



Color

Color option for feet: white,blue, pink,Tiffany blue ,yellow. MOQ500pcs for each color,IV stands can stack together to save space.




Comprehensive Hospital Design Solutions




Toruscare has a complete industry chain targeting the global market.



Toruscare provides overall healthcare solutions.



Toruscare provides services for full coverage of target audiences.



Design and Development Excellence

At Toruscare, our strength lies in our robust design and development capabilities. Our dedicated team of engineers and designers continually push the boundaries of innovation, ensuring that our products are at the forefront of medical technology. From conceptualization to prototyping, we meticulously craft each product to meet the highest standards of functionality, durability, and ease of use.



Uncompromising Quality Control

Toruscare is deeply committed to quality, implementing stringent control measures at every manufacturing stage. Our state-of-the-art facilities, equipped with advanced testing and inspection tools, ensure each product meets rigorous standards. This dedication to excellence has earned the trust and loyalty of healthcare providers worldwide.



Global Reach, Local Touch

As a global company, Toruscare understands the diverse needs of different markets. We pride ourselves on our ability to deliver customized solutions that cater to the unique requirements of healthcare providers across the world. Our extensive distribution network and strategic partnerships ensure timely delivery and exceptional after-sales support, making us a trusted partner in the global healthcare community.



Care for All

TORUSCARE

Toruscare Sp. z o.o

Adress: UL. SPACEROWA 1C, 41-253 CZELADŹ, POLAND

Email: export@toruscare.com

Phone: +351 910200400